

Product range PAX® & Artuvetrin®




Why consider allergy?

Today, one out of five visits to a veterinarian is about allergy and skin problems. Prevalence has increased rapidly over the past 10 years and is expected to continue to grow. Allergy is one of the world's most common chronic diseases affecting millions of dogs, cats and horses. It is estimated that up to 10% of the dogs and cats may suffer from allergic disease.

Pruritus is the most common symptom of allergy. Cutaneous lesions such as erythema, excoriations and self-induced alopecia also reflect the severity of the disease. Respiratory symptoms (allergic rhinitis, asthma), gastrointestinal symptoms (vomiting, diarrhea, gastritis) and ocular symptoms (conjunctivitis) may also appear. The diagnosis requires a systematic approach based on a full clinical examination, detailed patient history and a workup that excludes other potential causes.

The only causal treatment available is allergen-specific immunotherapy, a custom-made vaccine that treats the cause, not only the symptoms. It has proven efficacy (up to 70% in dogs) and safety, and avoids the need of long-term symptomatic medications with possible side effects.



Why serological testing?

Serological tests are a quick and easy way to identify the allergens that can be responsible for the allergic reaction.

Besides representing a quick and easy test, it has also the following advantages:

- Requires minimal patient preparation
- Result easy to interpret and with additional treatment advice
- No adverse reactions
- No need for referral
- Standardized procedure
- Excellent reproducibility
- Not influenced by existing skin pathologies

Following serological tests are available:

PAX® Complete

- Environmental
- Food
- Environmental + Food

Fcε-receptor Test

- Seasonal allergens
- Perennial allergens
- Insects / hymenoptera
- Feathers / hair / epithelia
- Flea saliva
- Mediterranean panel
- Malassezia

You can find more information about each one of them in the following pages.



In case of any question, please contact our dermatology team:

Tel: +49 971 7202 0 | Email: derma@laboklin.com

Why PAX[®] - Pet Allergy Xplorer?

Pet Allergy Xplorer (PAX) is the first commercial serological IgE-specific test that uses allergen extracts and molecular components to identify which allergens are affecting pets.

Molecular allergology is a state-of-the-art approach to the detection of sensitisations, whereby defined single allergen components are used for the determination of specific IgE in place of traditionally-used allergen extracts. The molecular components are purified or recombinant proteins that provide a higher level of standardisation than allergen extracts and enable a more precise identification of IgE sensitisations.

PAX advantages:

- 1 First quantitative multiplex macroarray specifically designed for companion animals
- 2 Over 200 allergens included = lower testing cost per allergen
- 3 Fully automated process = higher level of standardization (same result if tested multiple times)
- 4 With CCD blocking and 2 blocking efficiency detectors
- 5 Only 0.5 ml of serum needed
- 6 Expected increase in serological test sensitivity due to a higher concentration of molecular allergens
- 7 Identification of "primary" sensitizing allergens
- 8 Identification of allergen cross-reactivities
- 9 Selection of relevant allergens for specific immunotherapy

What are CCDs and why is important to block them?

CCDs are cross-reactive carbohydrate determinants – the carbohydrate chains found in glycoproteins. CCDs are part of the structure of many allergy-causing proteins, especially plant-based allergens like pollen. During an allergic reaction, IgE is produced against the carbohydrate chains as well as the allergen proteins.

Studies have confirmed that this is the case for 30% of humans, dogs and cats. Further research shows that specific IgEs against CCD are not clinically relevant. According to the latest studies, removing or blocking CCD specific IgE improves the correlation between in vitro tests and intradermal tests. Blocking CCDs simply means that the specificity of the in vitro test is enhanced.


PAX is the first IgE test to include not only CCD blocking, but also 2 blocking efficiency detectors that will evaluate and guarantee the completeness of the blocking.

EXTRACT




Allergen extract

MOLECULAR COMPONENTS




Der f 2
NPC2 family



Der f 3
trypsin



Der f 10
tropomyosin



Der f 15
chitinase



Der f 1
cysteine protease


PAX[®] - Pet Allergy Xplorer

Sample required:

0.5 ml
of serum

Available for:

Dog
Cat
Horse



Tested allergen components and extracts

	Common name	Scientific name	Extracts & Components
Grass Pollen	Bermuda grass	<i>Cynodon dactylon</i>	Cyn d * Cyn d 1
	Orchard grass	<i>Dactylis glomerata</i>	Dac g *
	Meadow fescue	<i>Festuca pratensis</i>	Fes p *
	Perennial ryegrass	<i>Lolium perenne</i>	Lol p 1
	Timothy	<i>Phleum pratense</i>	Phl p 1
			Phl p 2
			Phl p 5.0101
			Phl p 6
	Kentucky blue grass	<i>Poa pratensis</i>	Phl p 7
Phl p 12			
Poa p *			
Ryegrass, cultivated	<i>Secale cereale</i>	Sec c_pollen *	
Tree Pollen	Alder	<i>Alnus glutinosa</i>	Aln g *
			Aln g 1
			Aln g 4
	Silver birch	<i>Betula verrucosa</i>	Bet v *
			Bet v 1
			Bet v 2
			Bet v 6
Hazel	<i>Corylus avellana</i>	Cor a_pollen * Cor a 1.0103	

	Common name	Scientific name	Extracts & Components
Tree Pollen	Cypress	<i>Cupressus sempervirens</i>	Cup s *
	Beech	<i>Fagus sylvatica</i>	Fag s 1
	Ash	<i>Fraxinus excelsior</i>	Fra e *
			Fra e 1
	Privet	<i>Ligustrum vulgare</i>	Lig v *
	Olive tree	<i>Olea Europaea</i>	Ole e *
			Ole e 1
			Ole e 7
			Ole e 9
	London plane tree	<i>Platanus acerifolia</i>	Pla a 1
Pla a 2			
Pla a 3			
Cottonwood	<i>Populus nigra</i>	Pop n *	
Elm	<i>Ulmus campestris</i>	Ulm c *	
Weed Pollen	Ragweed	<i>Ambrosia artemisiifolia</i>	Amb a *
			Amb a 1
			Amb a 4
	Mugwort	<i>Artemisia vulgaris</i>	Art v *
			Art v 1.0101 Art v 3.0201
Lamb's quarter	<i>Chenopodium album</i>	Che a *	
		Che a 1	

	Common name	Scientific name	Extracts & Components	
Weed Pollen	Wall pellitory	<i>Parietaria judaica</i>	Par j *	
			Par j 2	
	Ribwort / Plantain	<i>Plantago lanceolata</i>	Pla l *	
			Pla l 1	
	Dock / Sorrel	<i>Rumex crispus / acetosella</i>	Rum c / * Rum a	
	Russian thistle	<i>Salsola kali</i>	Sal k *	
Sal k 1				
Nettle	<i>Urtica dioica</i>	Urt d *		
Danders & Epithelia	Cattle	<i>Bos domesticus</i>	Bos d 2	
	Dog	<i>Canis familiaris</i>	Can f 1	
			Can f 2	
			Can f 3	
			Can f 4	
			Can f 6	
				Can f_maleurine (including Can f 5) * Can f Fel d 1 like
	Guinea pig	<i>Cavia porcellus</i>	Cav p 1	
	Horse	<i>Equus caballus</i>	Equ c 1	
			Equ c 3 Equ c 4	

* Extract

	Common name	Scientific name	Extracts & Components	
Danders & Epithelia	Cat	<i>Felis domesticu</i>	Fel d 1	
			Fel d 2	
			Fel d 4	
			Fel d 7	
Mouse	<i>Mus musculus</i>	Mus m 1		
Rabbit	<i>Oryctolagus cuniculus</i>	Ory c 1		
		Ory c 2		
		Ory c 3		
Acarus siro	<i>Acarus siro</i>	Aca s *		
German cockroach	<i>Blattella germanica</i>	Bla g 1		
		Bla g 2		
		Bla g 4		
		Bla g 5		
		Bla g 9		
Flea	<i>Ctenocephalides felis</i>	Cte f 1		
Dermatophagoides pteronyssinus	<i>Dermatophagoides farinae</i>	Der f *		
		Der f 1		
		Der f 2		
		Der f 15		
Dermatophagoides pteronyssinus	<i>Dermatophagoides pteronyssinus</i>	Der p *		
		Der p 1		
		Der p 2		
		Der p 5		
Dermatophagoides pteronyssinus	<i>Dermatophagoides pteronyssinus</i>	Der p 7		
		Der p 10		
		Der p 11		
		Der p 20		
Dermatophagoides pteronyssinus	<i>Dermatophagoides pteronyssinus</i>	Der p 21		
		Der p 23		
		Glycyphagus domesticus	<i>Glycyphagus domesticus</i>	Gly d 2
		Lepidoglyphus destructor	<i>Lepidoglyphus destructor</i>	Lep d *
Lep d 2				
Tyrophagus putrescentiae	<i>Tyrophagus putrescentiae</i>	Tyr p *		
		Tyr p 2		
Moulds & Yeasts	Alternaria alternata	<i>Alternaria alternata</i>	Alt a *	
			Alt a 1	
			Alt a 6	
	Aspergillus fumigatus	<i>Aspergillus fumigatus</i>	Asp f *	
			Asp f 1	
			Asp f 3	
			Asp f 4	
			Asp f 6	
	Cladosporium herbarum	<i>Cladosporium herbarum</i>	Cla h *	
	Cladosporium herbarum	<i>Cladosporium herbarum</i>	Cla h 8	
	Malassezia pachydermatis	<i>Malassezia pachydermatis</i>	Mala p *	
	Malassezia sympodialis	<i>Malassezia sympodialis</i>	Mala s 1	
			Mala s 9	
			Mala s 5	
Mala s 6				
Mala s 11				

	Common name	Scientific name	Extracts & Components
Insect Venoms	Honey bee venom	<i>Apis mellifera</i>	Api m *
			Api m 1
			Api m 2
			Api m 3
			Api m 5
	Api m 10		
Long-headed wasp venom	<i>Dolichovespula spp.</i>	Dol spp *	
Paper wasp venom	<i>Polistes dominulus</i>	Pol d *	
		Pol d 5	
Fire ant venom	<i>Solenopsis richteri & Solenopsis invicta</i>	Sol spp *	
Common wasp venom	<i>Vespula vulgaris</i>	Ves v *	
		Ves v 1	
		Ves v 5	
Oat	<i>Avena sativa</i>	Ave s *	
Buckwheat	<i>Fagopyrum esculentum</i>	Fag e *	
		Fag e 2	
Sunflower seed	<i>Helianthus annuus</i>	Hel a *	
Barley	<i>Hordeum vulgare</i>	Hor v *	
Rice	<i>Oryza sativa</i>	Ory s	
		Saatroggen	
Millet	<i>Panicum miliaceum</i>	Pan m *	
Rye, cultivated	<i>Secale cereale</i>	Sec c_flour *	
Wheat	<i>Triticum aestivum</i>	Tri a *	
		Erdnuss	
		Soja	
		Linse	
Corn, cereal	<i>Zea mays</i>	Zea m *	
		Zea m 14	
		Zea m_GBSSI	
Apple	<i>Malus domestica</i>	Mal d 1	
		Mal d 2	
		Mal d 3	
Peanut	<i>Arachis hypogaea</i>	Ara h 1	
		Ara h 2	
		Ara h 3	
		Ara h 5	
		Ara h 6	
		Ara h 8	
Ara h 9			
Ara h 15			
Soy	<i>Glycine max</i>	Gly m *	
		Gly m 4	
		Gly m 5	
		Gly m 6	
		Gly m 8	
		Gly m 8	
Lentil	<i>Lens culinaris</i>	Len c *	
		Len c 1	
		Len c 2	
		Len c 3	
Pea	<i>Pisum sativum</i>	Pis s *	
		Pis s 1	
		Pis s 2	
		Pis s 3	

	Common name	Scientific name	Extracts & Components
Foods	Cow's milk	<i>Bos domesticus</i>	Bos d_milk *
			Bos d 4
			Bos d 5
			Bos d 8
	Egg white	<i>Gallus domesticus</i>	Gal d_white *
			Gal d 1
			Gal d 2
			Gal d 3
			Gal d 4
	Egg yolk	<i>Gallus domesticus</i>	Gal d_yolk *
			Gal d 5
	Beef	<i>Bos domesticus</i>	Bos d_meat *
			Bos d 6
			Bos d 7
			Bos d_ACTA1
			Bos d_LDHA
	Horse	<i>Equus caballus</i>	Equ c_meat *
	Rabbit	<i>Oryctolagus spp.</i>	Ory_meat *
			Ory c_CKM
			Ory c_GAPDH
			Ory c_PGM1
			Ory c_PKM
			Ory c_TPI1
	Lamb	<i>Ovis aries</i>	Ovi a_meat *
			Ovi a_IgG
	Pig	<i>Sus domesticus</i>	Sus d_meat *
Sus d 1			
Chicken	<i>Gallus domesticus</i>	Gal d_meat *	
		Gal d 7	
		Gal d 9	
		Gal d_PKM	
Turkey	<i>Meleagris gallopavo</i>	Mel g *	
Mealworm	<i>Tenebrio molitor</i>	Ten m *	
Herring, Atlantic	<i>Clupea harengus</i>	Clu h *	
		Clu h 1	
Cod, Atlantic	<i>Gadus morhua</i>	Gad m *	
		Gad m 1	
		Gad m 2+3	
		Gad m 4	
Salmon, Atlantic	<i>Salmo salar</i>	Sal s *	
		Sal s 1	
		Sal s 2	
		Sal s 3	
		Sal s 4	
		Sal s 6	
Sal s 7			
Sal s 8			
Mackerel, Atlantic	<i>Scomber scombrus</i>	Sco s *	
		Sco s 1	
Tuna	<i>Thunnus albacares</i>	Thu a *	
		Thu a 1	
Carrot	<i>Daucus carota</i>	Dau c *	
		Dau c 1	
Tomato	<i>Solanum lycopersicum</i>	Sola l *	
		Sola l 6	
Potato	<i>Solanum tuberosum</i>	Sol t *	
		Sol t 2	
		Sol t_GBSSI	

* Extract

PAX[®] Complete

- Environmental + Food
- Environmental
- Food

LABOKLIN
LABOR FÜR KLINISCHE DIAGNOSTIK GMBH & CO. KG
Südbödenstr. 4
97688 Bad Kissingen
Tel: +49 9717 20 20
info@laboklin.com

Owner Name: _____
Animal Name: _____
Species: _____
Test type: PAX Environmental
Date: _____
Lab Number: _____
Order Number: _____
Customer number: _____



PAX
Environmental

Summary on detectable sensitizations

<p>GRASS POLLENS</p> <p>Bermuda grass 31 Kentucky blue grass 22 Maize pollen 25 Orchard grass 22 Panicum polyanthemum 25 Rye, cultivated 25 Timothy grass 25</p> <p>WEED POLLENS</p> <p>Beet 21 English Plantain 26 Lamb's quarter 23 Nettle 23 Pillwort 23 Ragweed 28 Russet-top 23</p> <p>TREE POLLENS</p> <p>Alder 26 Ash 23 Birch 22 Beech 22 Cypress 22 Elaeagnus 24 Hornbeam 22 Olive tree 23 Pine 23 Pistachio 25 Spruce 22</p>	<p>MITES & COCKROACHES</p> <p>Acarus spp. 21 D. farinae 30 D. pteronyssinus 22 Glycyphagus domesticus 20 Lepidoglyphus destructor 29 Tyrophagus putrescentiae 29 German cockroach 28 American cockroach 18 Flea 18</p> <p>MOULDS & YEASTS</p> <p>Aspergillus fumigatus 25 Aspergillus terreus 24 Cladosporium herbarum 24 Penicillium 63</p> <p>EPITHELIA</p> <p>Cat 24 Cattle 19 Guinea pig 21 Horse 21 Muscle 20 Rabbit 22</p> <p>INSECT VENOMS</p> <p>Fire ant venom 25 Honey bee venom 21 Long-headed wasp venom 21 Paper wasp venom 22</p>
---	---

The results above represent a summary of the sensitizations detected with extracts and components for each allergen. Detailed results can be found in the following pages.

Highest measured IgE concentration per allergen group

< 28.00 ng/mL	28.00-99.99 ng/mL	100.00-399.99 ng/mL	400.00-799.99 ng/mL	≥ 800.00 ng/mL
class 0	class 1	class 2	class 3	class 4

Class 0 - Values below 28 ng/ml are usually not considered relevant for immunotherapy. However, do note that the serum level of allergen-specific IgE is only one of the factors of IgE's potential to cause clinical signs of allergy. It is likely that, in some cases, allergies might occur with levels of specific IgE below that cutoff value. If allergy signs are associated with low levels of IgE (< 28 ng/ml), equine allergens suspected to be relevant for the patient, then immunotherapy could be considered.

LABOKLIN
LABOR FÜR KLINISCHE DIAGNOSTIK GMBH & CO. KG

Owner Name: _____
Animal Name: _____
Species: _____
Test type: PAX Environmental
Date: _____
Lab Number: _____
Order Number: _____
Customer number: _____



PAX
Environmental

Summary on detectable sensitizations

<p>GRASS POLLENS</p> <p>Bermuda grass Cynodon dactylon Cyn 4 18 Kentucky blue grass Poa pratensis Psa 9 22 Maize pollen Festuca pratensis Fes 9 22 Orchard grass Dactylis glomerata Dac 9 22 Panicum polyanthemum Lolium perenne Lol p 1 25 Rye, cultivated Secale cereale Sec C Lyden 25 Timothy grass Phleum pratense Pfl p 1 22</p> <p>WEED POLLENS</p> <p>Beet / Sorrel Rumex acetosella / crispus Rum a, Rum c 22 English Plantain Plantago lanceolata Pla 1 26 Lamb's quarter Chenopodium album Cha 4 24 Ragweed Artemisia vulgaris Art v 19 Russet-top Ambrosia artemisiifolia Art a 1 22 Nettle Urtica dioica Ur t 6 23 Pillwort Parietaria judicaria Par j 21 Ragweed Ambrosia artemisiifolia Art a 1 22 Russet-top Salvia sal Sal 1 23</p>	<p>MITES & COCKROACHES</p> <p>Acarus spp. 21 D. farinae 30 D. pteronyssinus 22 Glycyphagus domesticus 20 Lepidoglyphus destructor 29 Tyrophagus putrescentiae 29 German cockroach 28 American cockroach 18 Flea 18</p> <p>MOULDS & YEASTS</p> <p>Aspergillus fumigatus 25 Aspergillus terreus 24 Cladosporium herbarum 24 Penicillium 63</p> <p>EPITHELIA</p> <p>Cat 24 Cattle 19 Guinea pig 21 Horse 21 Muscle 20 Rabbit 22</p> <p>INSECT VENOMS</p> <p>Fire ant venom 25 Honey bee venom 21 Long-headed wasp venom 21 Paper wasp venom 22</p>
--	---

The results above represent a summary of the sensitizations detected with extracts and components for each allergen. Detailed results can be found in the following pages.

Highest measured IgE concentration per allergen group

< 28.00 ng/mL	28.00-99.99 ng/mL	100.00-399.99 ng/mL	400.00-799.99 ng/mL	≥ 800.00 ng/mL
class 0	class 1	class 2	class 3	class 4

Class 0 - Values below 28 ng/ml are usually not considered relevant for immunotherapy. However, do note that the serum level of allergen-specific IgE is only one of the factors of IgE's potential to cause clinical signs of allergy. It is likely that, in some cases, allergies might occur with levels of specific IgE below that cutoff value. If allergy signs are associated with low levels of IgE (< 28 ng/ml), equine allergens suspected to be relevant for the patient, then immunotherapy could be considered.

LABOKLIN
LABOR FÜR KLINISCHE DIAGNOSTIK GMBH & CO. KG

Owner Name: _____
Animal Name: _____
Species: _____
Test type: PAX Environmental
Date: _____
Lab Number: _____
Order Number: _____
Customer number: _____



PAX
Environmental

Interpretation - Support

Pellitory

- This patient has a sensitization to pellitory pollen.
- Associated allergic signs are generally worse during the pellitory pollen season, which varies from May to October depending upon the geographical location.
- Allergen-specific immunotherapy is recommended for pellitory pollen sensitization, if the corresponding clinical signs occur.
- Pari j 2 is an allergen from the wall pellitory (Parietaria judicaria). It is a member of the mUP2 allergen family.
- The potential for cross-reactions with most other allergens of this family is considered low.
- In humans, Pari j 2 is a highly specific marker for pellitory pollen sensitization; at this time, it is not known if the same occurs in animals.

D. farinae

- This patient has a sensitization to house dust mites.
- Associated allergic signs are generally year-round, but house dust mites are known to proliferate during times of high humidity and temperature.
- There is a known cross-reactivity between allergens of house dust and storage mite species, as well as between those of Dermatophagoides farinae and D. pteronyssinus.
- Allergen-specific immunotherapy is recommended for house dust mite sensitization, if the corresponding clinical signs occur.

Tyrophagus putrescentiae

- This patient has a sensitization to storage mites.
- Associated allergic signs are generally year-round, but storage mites are known to proliferate during times of high humidity and temperature.
- There is a known cross-reactivity between allergens of house dust and storage mite species.
- Allergen-specific immunotherapy is recommended for storage mite sensitization, if the corresponding clinical signs occur.

American cockroach

- This patient has a sensitization to American cockroach.
- Associated allergic signs are generally year-round.
- There is a known cross-reactivity between allergens of different cockroach species.
- Allergen-specific immunotherapy is recommended for cockroach sensitization, if the corresponding clinical signs occur.
- Pari 7 is an allergen of the American cockroach (Parietaria americana). It belongs to the family of cockroach Group 7 allergens (Pteromyzins).
- Pari 7 is a major allergen of humans sensitized to this cockroach; at this time it is not known if this is also the case in animals.
- The potential for cross-reactions of Pari 7 with other tropomyosins present in insects, nematodes, and ingested seafood is very high.
- Sensitization to tropomyosin likely represents a cross-reactivity to nematodes, such as Tricostema. At this time, the benefit of immunotherapy for patient sensitized to cross-reactive tropomyosins has not been studied. However, in the absence of additional sensitizations, immunotherapy to tropomyosin-containing house dust mites could be considered.

DISCLAIMER: THIS INTERPRETATION GUIDANCE IS A TOOL TO SUPPORT VETERINARIANS IN THE INTERPRETATION OF PAX RESULTS. PAX RECOMMENDATIONS DO NOT REPLACE THE JUDGMENT OF A VETERINARIAN.

LABOKLIN
LABOR FÜR KLINISCHE DIAGNOSTIK GMBH & CO. KG

Owner Name: _____
Animal Name: _____
Species: _____
Test type: PAX Environmental, Insects + Food
Date: _____
Lab Number: _____
Order Number: _____
Customer number: _____



PAX
Environmental

Summary on detectable sensitizations

<p>GRASS POLLENS</p> <p>Bermuda grass 20 Kentucky blue grass 23 Maize pollen 25 Orchard grass 22 Panicum polyanthemum 25 Rye, cultivated 25 Timothy grass 25</p> <p>WEED POLLENS</p> <p>Beet 21 English Plantain 26 Lamb's quarter 23 Nettle 23 Pillwort 23 Ragweed 28 Russet-top 23</p> <p>TREE POLLENS</p> <p>Alder 26 Ash 23 Birch 22 Beech 22 Cypress 22 Elaeagnus 24 Hornbeam 22 Olive tree 23 Pine 23 Pistachio 25 Spruce 22</p>	<p>MITES & COCKROACHES</p> <p>Acarus spp. 21 D. farinae 30 D. pteronyssinus 22 Glycyphagus domesticus 20 Lepidoglyphus destructor 29 Tyrophagus putrescentiae 29 German cockroach 28 American cockroach 18 Flea 18</p> <p>MOULDS & YEASTS</p> <p>Aspergillus fumigatus 25 Aspergillus terreus 24 Cladosporium herbarum 24 Penicillium 63</p> <p>EPITHELIA</p> <p>Cat 24 Cattle 19 Guinea pig 21 Horse 21 Muscle 20 Rabbit 22</p> <p>BITING INSECTS</p> <p>Culicoides 23 Deer fly 25 House fly 22 Horse fly 25 Stable fly 25</p> <p>INSECT VENOMS</p> <p>Fire ant venom 25 Honey bee venom 21 Long-headed wasp venom 20 Paper wasp venom 22</p> <p>OTHER</p> <p>Lactose 25</p>
---	---

The results above represent a summary of the sensitizations detected with extracts and components for each allergen. Detailed results can be found in the following pages.

Highest measured IgE concentration per allergen group

< 28.00 ng/mL	28.00-99.99 ng/mL	100.00-399.99 ng/mL	400.00-799.99 ng/mL	≥ 800.00 ng/mL
class 0	class 1	class 2	class 3	class 4

Class 0 - Values below 28 ng/ml are usually not considered relevant for immunotherapy. However, do note that the serum level of allergen-specific IgE is only one of the factors of IgE's potential to cause clinical signs of allergy. It is likely that, in some cases, allergies might occur with levels of specific IgE below that cutoff value. If allergy signs are associated with low levels of IgE (< 28 ng/ml), equine allergens suspected to be relevant for the patient, then immunotherapy could be considered.

LABOKLIN
LABOR FÜR KLINISCHE DIAGNOSTIK GMBH & CO. KG

Owner Name: _____
Animal Name: _____
Species: _____
Test type: PAX Environmental, Insects + Food
Date: _____
Lab Number: _____
Order Number: _____
Customer number: _____



PAX
Environmental

Summary on detectable sensitizations

<p>GRASS POLLENS</p> <p>Bermuda grass Cynodon dactylon Cyn 4 18 Kentucky blue grass Poa pratensis Psa 9 22 Maize pollen Festuca pratensis Fes 9 22 Orchard grass Dactylis glomerata Dac 9 22 Panicum polyanthemum Lolium perenne Lol p 1 25 Rye, cultivated Secale cereale Sec C Lyden 25 Timothy grass Phleum pratense Pfl p 1 22</p> <p>WEED POLLENS</p> <p>Beet / Sorrel Rumex acetosella / crispus Rum a, Rum c 24 English Plantain Plantago lanceolata Pla 1 26 Lamb's quarter Chenopodium album Cha 4 24 Ragweed Artemisia vulgaris Art v 19 Russet-top Ambrosia artemisiifolia Art a 1 22 Nettle Urtica dioica Ur t 6 23 Pillwort Parietaria judicaria Par j 21 Ragweed Ambrosia artemisiifolia Art a 1 22 Russet-top Salvia sal Sal 1 23</p>	<p>MITES & COCKROACHES</p> <p>Acarus spp. 21 D. farinae 30 D. pteronyssinus 22 Glycyphagus domesticus 20 Lepidoglyphus destructor 29 Tyrophagus putrescentiae 29 German cockroach 28 American cockroach 18 Flea 18</p> <p>MOULDS & YEASTS</p> <p>Aspergillus fumigatus 25 Aspergillus terreus 24 Cladosporium herbarum 24 Penicillium 63</p> <p>EPITHELIA</p> <p>Cat 24 Cattle 19 Guinea pig 21 Horse 21 Muscle 20 Rabbit 22</p> <p>BITING INSECTS</p> <p>Culicoides 23 Deer fly 25 House fly 22 Horse fly 25 Stable fly 25</p> <p>INSECT VENOMS</p> <p>Fire ant venom 25 Honey bee venom 21 Long-headed wasp venom 20 Paper wasp venom 22</p> <p>OTHER</p> <p>Lactose 25</p>
--	---

The results above represent a summary of the sensitizations detected with extracts and components for each allergen. Detailed results can be found in the following pages.

Highest measured IgE concentration per allergen group

< 28.00 ng/mL	28.00-99.99 ng/mL	100.00-399.99 ng/mL	400.00-799.99 ng/mL	≥ 800.00 ng/mL
class 0	class 1	class 2	class 3	class 4

Class 0 - Values below 28 ng/ml are usually not considered relevant for immunotherapy. However, do note that the serum level of allergen-specific IgE is only one of the factors of IgE's potential to cause clinical signs of allergy. It is likely that, in some cases, allergies might occur with levels of specific IgE below that cutoff value. If allergy signs are associated with low levels of IgE (< 28 ng/ml), equine allergens suspected to be relevant for the patient, then immunotherapy could be considered.

LABOKLIN
LABOR FÜR KLINISCHE DIAGNOSTIK GMBH & CO. KG

Owner Name: _____
Animal Name: _____
Species: _____
Test type: PAX Environmental, Insects + Food
Date: _____
Lab Number: _____
Order Number: _____
Customer number: _____



PAX
Environmental

Interpretation - Support

D. pteronyssinus

- This patient has a sensitization to house dust mites.
- Associated allergic signs are generally year-round, but house dust mites are known to proliferate during times of high humidity and temperature.
- There is a known cross-reactivity between allergens of house dust and storage mite species.
- Allergen-specific immunotherapy is recommended for house dust mite sensitization, if the corresponding clinical signs occur.
- Der p 1 is an allergen from the Dermatophagoides pteronyssinus house dust mite. It is a member of the mite Group 1 allergen family (lipid transfer proteins).
- Der p 1 is a major allergen of humans sensitized to this house dust mite; it has also been shown to be an allergen in mite-sensitized dogs. At this time we do not know if this is the case in cats and horses.
- The potential for cross-reactions with other Group 1 mite allergens is variable, from very high (Der f 1 and Der m 1) to low (Ty r p 1, Bix 1, Aca 1, Sar 1).

Mosquito

- This patient has a sensitization to mosquito.
- Clinical signs are due to a combination of immediate and late-phase reactions that follow the bites of insects; signs are generally worse in the summer and fall.
- There is a high likelihood of cross-reactivity between salivary allergens of different biting insect species, especially those that are closely related.
- The avoidance of insect bites via multifaceted insect control measures is the most effective treatment option for insect bite hypersensitivity. Immunotherapy is recommended for these insects, if the corresponding clinical signs occur.

Honey bee venom

- This patient has a sensitization to honey bee venom.
- Clinical signs follow stings from honey bees, and involve a combination of immediate and late-phase reactions at the site of stings; systemic signs also can occur.
- There is a known cross-reactivity between allergens of the honeybee and wasp venoms. In humans, cross-reactivity has been demonstrated between the bee venoms and vespid venoms.
- Allergen-specific immunotherapy is currently not available for honey bee venom sensitization; the treatment is symptomatic.
- Api m 1 is an allergen from the honey bee (Apis mellifera). It is the honey bee venom's phospholipase A2.
- Api m 1 is a CD2-carrying major allergen of humans sensitized to honey bee venom; our preliminary results suggest that most dogs sensitized to the honey bee venom also are sensitized to Api m 1. At this time we do not know if this is the case in cats and horses.
- Api m 1 is considered a marker for sensitization to honey bee and bumble bee; phospholipase A2 are not cross-reactive with phospholipase A1 from vespid (e.g., Vex s 1).

DISCLAIMER: THIS INTERPRETATION GUIDANCE IS A TOOL TO SUPPORT VETERINARIANS IN THE INTERPRETATION OF PAX RESULTS. PAX RECOMMENDATIONS DO NOT REPLACE THE JUDGMENT OF A VETERINARIAN.

The PAX[®] Complete panel specifies the individual result for each tested allergen extract and component. It includes a detailed recommendation including possible cross-reactivities and a support interpretation for each positive tested allergen.

Artuvetrin® Skin Test

The Artuvetrin® Skin Test allows you to determine to which substances a dog is allergic, by injecting small amounts of different allergens in the skin (intradermal test).



Latin	English	Mixtures
-------	---------	----------

Pollen of grasses		
-------------------	--	--

- | | | |
|--------------------------|----------------------|---|
| 1. Anthoxanthum odoratum | Vernal grass, sweet | Grass pollen mixture: <ul style="list-style-type: none"> Sweet vernal grass, bermuda grass, orchard grass, velvet grass and timothy grass |
| 2. Cynodon dactylon | Bermuda grass | |
| 3. Dactylis glomerata | Orchard grass | |
| 4. Festuca pratensis | Fescue, meadow | |
| 5. Holcus lanatus | Velvet grass | |
| 6. Lolium perenne | Rye grass, perennial | |
| 7. Phleum pratense | Timothy grass | |
| 8. Poa pratensis | Blue grass, Kentucky | |

Pollen of trees and hedges		
----------------------------	--	--

- | | | |
|---------------------------|-------------------|---|
| 1. Acer negundo | Box elder | Tree pollen mixture 1: <ul style="list-style-type: none"> Alder, birch and hazel Tree pollen mixture 2: <ul style="list-style-type: none"> Beech, oak and elm |
| 2. Corylus avellana | Hazel | |
| 3. Cupressus sempervirens | Cypress, Italian | |
| 4. Fagus sylvatica | Beech, European | |
| 5. Olea europaea | Olive | |
| 6. Platanus occidentalis | Sycamore, eastern | |
| 7. Populus alba | Poplar, white | |

Pollen of weeds and flowers		
-----------------------------	--	--

- | | | |
|---------------------------|--------------------------|--|
| 1. Ambrosia elatior | Ragweed, common | Weed pollen mixture 1: <ul style="list-style-type: none"> Mugwort, English plantain, dandelion and nettle Weed pollen mixture 2: <ul style="list-style-type: none"> Mugwort and nettle Weed pollen mixture 3: <ul style="list-style-type: none"> Lamb's quarters, English plantain, red or sheep sorrel, goldenrod and dandelion |
| 2. Artemisia vulgaris | Mugwort, common | |
| 3. Chenopodium album | Lamb's quarters | |
| 4. Parietaria officinalis | Pellitory, wall or erect | |
| 5. Plantago lanceolata | Plantain, English | |
| 6. Rumex acetosella | Sorrel, red or sheep | |
| 7. Salsola kali | Russian Thistle | |
| 8. Solidago virgaurea | Goldenrod | |

Latin	English	Mixtures
-------	---------	----------

Pollen of crops		
-----------------	--	--

- | | | |
|-------------------|-------------|-----|
| 1. Brassica napus | Rape | N/A |
| 2. Secale cereale | Rye, common | |

Mites		
-------	--	--

- | | | |
|-----------------------------------|-----------------|-----|
| 1. Acarus siro | Grain mite | N/A |
| 2. Dermatophagoides farinae | House dust mite | |
| 3. Dermatophagoides pteronyssinus | House dust mite | |
| 4. Lepidoglyphus destructor | Hay mite | |
| 5. Tyrophagus putrescentiae | Copra mite | |

Insects		
---------	--	--

- | | | |
|--------------------------|---------------|-----|
| 1. Aedes communis | Mosquito | N/A |
| 2. Culex | Culex | |
| 3. Culicoides | Biting Midges | |
| 4. Musca domestica | House fly | |
| 5. Periplaneta americana | Cockroach | |

Epithelia		
-----------	--	--

- | | | |
|----------------------------|----------|--|
| 1. Canis familiaris | Dog | Epithelia mixture:
<ul style="list-style-type: none"> Duck, goose and chicken |
| 2. Felis domesticus | Cat | |
| 3. Melopsittacus undulatus | Parakeet | |
| 4. Mus musculus | Mouse | |
| 5. Oryctolagus cuniculus | Rabbit | |
| 6. Ovis arie | Sheep | |

Fungi		
-------	--	--

- | | | |
|--------------------------|--|--|
| 1. Alternaria alternata | | Fungi mixture:
<ul style="list-style-type: none"> Alternaria alternata, Aspergillus fumigatus and Cladosporium herbarum |
| 2. Aspergillus fumigatus | | |
| 3. Cladosporium herbarum | | |

Yeast		
-------	--	--

- | | | |
|---------------------|--|-----|
| 1. Candida albicans | | N/A |
| 2. Malassezia | | |

This list is subjective to change. If you have questions about the availability, please contact derma@laboklin.com

Artuvetrin® allergen-specific immunotherapy (1 – 8 allergens)

Artuvetrin® Therapy is a licensed allergen-specific immunotherapy in Europe. It is a tailor-made treatment produced individually for each animal. Artuvetrin® Therapy has proven efficacy (75%) and safe. Therefore, it can be given for long periods. Each vial lasts for 10 months, and up to 8 allergens can be added to one vial.



Latin	English	Mixtures
Pollen of grasses		
<ol style="list-style-type: none"> 1. Anthoxanthum odoratum 2. Cynodon dactylon 3. Dactylis glomerata 4. Festuca pratensis 5. Holcus lanatus 6. Lolium perenne 7. Phleum pratense 8. Poa pratensis 	<p>Vernal grass, sweet Bermuda grass Orchard grass Fescue, meadow Velvet grass Rye grass, perennial Timothy grass Blue grass, Kentucky</p>	<p>Grass pollen mixture:</p> <ul style="list-style-type: none"> • Sweet vernal grass, bermuda grass, orchard grass, velvet grass and timothy grass
Pollen of trees and hedges		
<ol style="list-style-type: none"> 1. Acer negundo 2. Alnus glutinosa 3. Betula pendula 4. Corylus avellana 5. Cupressus sempervirens 6. Fagus sylvatica 7. Fraxinus excelsior 8. Ligustrum vulgare 9. Olea europaea 10. Platanus occidentalis 11. Populus alba 12. Salix viminalis 13. Ulmus americana 	<p>Maple, ash leaved Alder, black Birch Hazel Cypress, Italian Beech, European Ash European privet Olive Sycamore, eastern Poplar, white Willow Elm, American</p>	<p>Tree pollen mixture 1:</p> <ul style="list-style-type: none"> • Alder, birch and hazel <p>Tree pollen mixture 2:</p> <ul style="list-style-type: none"> • Beech, oak and elm <p>Tree pollen mixture 3:</p> <ul style="list-style-type: none"> • Ash-leaved maple, ash, white poplar and willow
Pollen of weeds and flowers		
<ol style="list-style-type: none"> 1. Ambrosia elatior 2. Artemisia vulgaris 3. Chenopodium album 4. Parietaria officinalis 5. Plantago lanceolata 6. Rumex acetosella 7. Salsola kali 8. Solidago virgaurea 9. Taraxacum officinale 10. Urtica dioica 	<p>Ragweed, common Mugwort, common Lamb's quarters Pellitory, wall or erect Plantain, English Sorrel, red or sheep Russian Thistle Goldenrod Dandelion, common Nettle, stinging</p>	<p>Weed pollen mixture 1:</p> <ul style="list-style-type: none"> • Mugwort, English plantain, dandelion and nettle <p>Weed pollen mixture 2:</p> <ul style="list-style-type: none"> • Mugwort and nettle <p>Weed pollen mixture 3:</p> <ul style="list-style-type: none"> • Lambs quarter, English plantain, red or sheep sorrel, goldenrod and dandelion

Latin	English	Mixtures
-------	---------	----------

Pollen of crops		
-----------------	--	--

- | | | |
|-------------------|-------------|-----|
| 1. Brassica napus | Rape | N/A |
| 2. Secale cereale | Rye, common | |

Mites		
-------	--	--

- | | | |
|-----------------------------------|-----------------|-----|
| 1. Acarus siro | Grain mite | |
| 2. Dermatophagoides farinae | House dust mite | |
| 3. Dermatophagoides pteronyssinus | House dust mite | N/A |
| 4. Lepidoglyphus destructor | Hay mite | |
| 5. Tyrophagus putrescentiae | Copra mite | |

Insects		
---------	--	--

- | | | |
|--------------------------|---------------|-----|
| 1. Aedes communis | Mosquito | |
| 2. Culex | Culex | |
| 3. Culicoides | Biting midges | N/A |
| 4. Musca domestica | House fly | |
| 5. Periplaneta americana | Cockroach | |

Epithelia		
-----------	--	--

- | | | |
|----------------------------|------------|--|
| 1. Canis familiaris | Dog | |
| 2. Cavia porcellus | Guinea pig | |
| 3. Felis domesticus | Cat | |
| 4. Melopsittacus undulatus | Parakeet | Epithelia mixture:
▪ Duck, goose and chicken |
| 5. Mus musculus | Mouse | |
| 6. Oryctolagus cuniculus | Rabbit | |
| 7. Ovis arie | Sheep | |

Fungi		
-------	--	--

- | | | |
|--------------------------|--|--|
| 1. Alternaria alternata | | Fungi mixture:
▪ Alternaria alternata,
Aspergillus fumigatus and
Cladosporium herbarum |
| 2. Aspergillus fumigatus | | |
| 3. Cladosporium herbarum | | |

Yeast		
-------	--	--

- | | | |
|---------------------|--|-----|
| 1. Candida albicans | | N/A |
| 2. Malassezia | | |

This list is subjective to change. If you have questions about the availability, please contact derma@laboklin.com

LABOKLIN

LABOKlin Labor für klinische Diagnostik GmbH & Co. KG

(Laboklin Laboratory for Clinical Diagnostics, a German limited partnership formed with a limited liability company)

D

Phone
Email
Internet

Steubenstr. 4
DE-97688 Bad Kissingen
Germany
+49 971 7 20 20
info@laboklin.com
www.laboklin.com