


Dermatology LABOKLIN

Sampling



Impression smear
Especially for oozing skin conditions




Skin scraping
Especially for the detection of parasites



Fine-needle biopsy
With or without aspiration
(pustules or masses)

Preparation



Smear



Air-dry



Stain
Rapid stain, e.g. Diff-Quick

Clinical and Microscopic Diagnosis



Pyoderma
Degenerate neutrophils and
bacterial cocci







Malassezia dermatitis
Malassezia and keratin flakes



Eosinophilic granuloma, cat
Eosinophils and squamous
epithelial cells




Sarcoptic mange



Pemphigus foliaceus
Acantholytic cells and neutrophils



Demodicosis
Demodex mites



Dermatophytosis
Fungal filaments



Cutaneous leishmaniasis
Macrophage with *Leishmania*

www.laboklin.com

As of: November 2019

D-97688 Bad Kissingen • Steubenstraße 4
Tel. +49-971/7 20 20 • Fax +49-971/6 85 46
E-mail: info@laboklin.com • www.laboklin.com